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




a xylem brand

EN

Thank you for purchasing this OPTi<sup>®</sup> Multiple Scale Digital Handheld Refractometer. In order to ensure that this product provides many years of service, please follow the guidance in this document.

In order to register your instrument with the manufacturer and print a warranty certificate, please visit:  
[www.bellinghamandstanley.com](http://www.bellinghamandstanley.com)

User Guide Code: 38-421'01PA-EN

## Unpacking the instrument

Check that all parts listed below are present and that no transit damage has occurred.

### Contents list

- 1 OPTi Digital Handheld Refractometer
- 2 AAA alkaline battery (LR03)
- 1 User guide

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## ! Safety precautions

### WARNINGS

Always check the relevant Material Safety Data Sheet for a sample before applying it to the refractometer. Wear appropriate protective equipment (PPE) when applying samples that may be harmful to the skin or eyes. Avoid unnecessary contamination of the refractometer by confining samples to the prism dish.

### CAUTION

This refractometer is a precision optical instrument and should be handled with care. Do not drop or subject the instrument to sharp knocks. The instrument housing and display panel areas are constructed from plastic materials that may suffer damage if contacted with aggressive solvents. For example, avoid contact with solvents such as acetone and certain aromatics.

Maintain your refractometer in a clean condition and avoid use and storage of the instrument outside the specified temperature range. Avoid dusty and high humidity environments and prolonged exposure to direct sunlight. Use the case provided to protect the instrument. Deterioration/loss of the display may be indicative of low battery power or low ambient temperature. Do not persist in using the instrument with low battery power. Check/replace the batteries as necessary.

### INTENDED USE

This product is for general laboratory, manufacturing and research use only and is not intended for any animal or human therapeutic or diagnostic use.

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All materials used to verify the performance of this OPTi refractometer were calibrated by Bellingham & Stanley Ltd., UKAS accredited calibration laboratory no. 0834, accredited to ISO/IEC 17025:2017.

1. Refractive index values for Water are obtained from "Revised Formulation for the Refractive Index of Water and Steam as a Function of Wavelength, Temperature and Density", adopted by the International Association for the Properties of Water and Steam (IAPWS) and available as part of NIST Standard Reference Database 10. Refractive indices calculated from the formulation are absolute refractive indices; conversion to refractive index against air requires division by the respective absolute refractive index of air (NIST Engineering Metrology Toolbox).

2. AG calibration fluids are prepared by mixing water with a soluble chemical of higher refractive index. The Brix value of the fluid is determined using a high accuracy digital refractometer at 20.0°C, which has been calibrated with a NIST traceable sucrose calibration standard of Brix value equal to the target value for the AG fluid according to a Laboratory Procedure QL-103. The AG fluid composition is finely adjusted to give a measured value to within ± 0.01 Brix (± 0.00002 RI) at 20.0°C.

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索引	應用	刻度	單位	範圍	解析度	準確度	溫度補償
01	基本	°Brix (ATC) (糖度含溫度補償)		0-95	0.1	±0.2	°Brix
02	基本	°Brix (糖度)		0-95	0.1	±0.2	None
03	基本	Refractive Index (ATC) (折射率含溫度補償)		1.33-1.53	0.0001	±0.0003	°Brix
04	基本	Refractive Index (折射率)		1.33-1.53	0.0001	±0.0003	None
05	汽車	AdBlue®/DEF (NOx reduction)	% Weight / Weight	0-40	0.1	±0.2	AUS32
06	汽車	Ethylene Glycol (乙二醇)	°C Freezing Point	0 to -50	1	±1	EG
07	汽車	Ethylene Glycol (乙二醇)	°F Freezing Point	30 to -40	1	±1	EG
08	汽車	Propylene Glycol (丙二醇)	°C Freezing Point	0 to -50	1	±1	PG
09	汽車	Propylene Glycol (丙二醇)	°F Freezing Point	30 to -40	1	±1	PG
10	食物飲料	°Butyro		0-100	0.1	±0.5	Butyro
11	食物飲料	Fructose (果糖)	% Weight / Weight	0-85	0.1	±0.2	°Brix
12	食物飲料	Glucose (葡萄糖)	% Weight / Weight	0-85	0.1	±0.2	°Brix
13	食物飲料	42 HFCS (高果糖玉米糖漿)	% Weight / Weight	0-95	0.1	±0.2	°Brix
14	食物飲料	55 HFCS (高果糖玉米糖漿)	% Weight / Weight	0-95	0.1	±0.2	°Brix
15	食物飲料	90 HFCS (高果糖玉米糖漿)	% Weight / Weight	0-95	0.1	±0.2	°Brix
16	食物飲料	Invert Sugar (轉化糖)	% Weight / Weight	0-85	0.1	±0.2	°Brix
17	食物飲料	Maltose (麥芽糖)	% Weight / Weight	0-60	0.1	±0.2	°Brix
18	食物飲料	Salinity (NaCl) (鹽度)	% Weight / Volume	0-28	0.1	±0.2	NaCl
19	食物飲料	Total Solids of Waste Milk (廢乳品總固形物)	%	5-15	0.1	±0.5	°Brix
20	食物飲料	Water in Honey (蜂蜜水份)	% Weight / Weight	10-30	0.1	±0.2	Honey
21	工業	Arbitrary (隨意)		0-95	0.1	±0.2	
22	工業	Calcium Chloride (氯化鈣)	% Weight / Weight	0-40	0.1	±0.2	CaCl <sub>2</sub>
23	工業	Ethanol (乙醇)	% Volume / Volume	0-20	0.5	±1	Ethanol
24	工業	Ethylene Glycol (乙二醇)	% Volume / Volume	0-60	0.1	±0.4	EG
25	工業	Ethylene Glycol (乙二醇)	% Weight / Weight	0-60	0.1	±0.4	EG
26	工業	FSII DiEGME (ASTM D 5006) (防凍劑)	% Volume / Volume	0.0-0.25	0.01	±0.02	°Brix
27	工業	Hydrogen Peroxide (過氧化氫)	% Weight / Weight	0-50	0.2	±0.5	°Brix
28	工業	Methanol (甲醇)	% Weight / Weight	0-40	1	±0.2	Meth
29	工業	Propylene Glycol (丙二醇)	% Volume / Volume	0-60	0.1	±0.4	PG
30	工業	Sodium Sulphate (硫酸鈉)	% Weight / Weight	0-22	0.1	±0.2	Na <sub>2</sub> SO <sub>4</sub>
31	工業	Starch (澱粉)	%	0-30	0.1	±0.2	°Brix
32	工業	Sulphuric Acid (Battery Acid) (硫酸)	Specific Gravity (d20/20)	1.000-1.501	0.001	±0.003	SA
33	工業	Urea (CRC data) (尿素)	% Weight / Weight	0-40	0.1	±0.2	AUS32
34	生命科學	Colostrum Quality (初乳品質)		Poor / PASS	Poor / PASS	±0.2	°Brix
35	生命科學	Seawater (Practical Salt Units) (海水)	Part Per Thousand	0-180	1	±1	NaCl
36	生命科學	Seawater (Practical Salt Units) (海水)	Specific Gravity (d20/20)	1.000-1.090	0.0005	±0.001	NaCl
37	生命科學	Serum Protein (血清蛋白)	g/100ml	0-30	0.1	±0.2	°Brix
38	生命科學	Urine (SG) Human (人類尿液)	Specific Gravity (d20/20)	1.000-1.050	0.0005	±0.0010	°Brix
39	生命科學	Urine (SG) Large Mammal (大型哺乳類尿液)	Specific Gravity (d20/20)	1.000-1.050	0.0001	±0.0010	°Brix
40	生命科學	Urine (SG) Small Mammal (小型哺乳類尿液)	Specific Gravity (d20/20)	1.000-1.050	0.0005	±0.0010	°Brix
41	酒類	°Baumé (波美比重)		0-50	0.1	±0.2	°Brix
42	酒類	°Zeiss (ABV)	% Volume / Volume	10-135	0.1	±0.5	°Brix
43	酒類	Alcohol Probable (AP)		0-22	0.1	±0.2	°Brix
44	酒類	KMW (Babo) (法國糖度)		0-25	1	±1	°Brix
45	酒類	Oechsle (German) (德國糖度)		30-130	1	±1	°Brix
46	酒類	Oechsle (Swiss) (瑞士糖度)		0-130	1	±1	°Brix
47	酒類	°Plato (柏拉圖度)		0-30	0.1	±0.2	°Brix
48	酒類	Mass Sugar (°Brix) (ATC)	% Weight / Weight	0-95	0.1	±0.2	°Brix
49	酒類	Mass Sugar (°Brix)	% Weight / Weight	0-95	0.1	±0.2	None
50	酒類	Wort (Sucrose Equivalent) (麥芽汁蔗糖當量)	Specific Gravity (d20/20)	1.000-1.120	0.0005	±0.001	°Brix

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This symbol is an internationally agreed indicator that the product bearing it should not be disposed of as general waste or garbage which might end up in landfill sites, but should instead be sent for special processing and/or recycling in those countries where appropriate legislation and facilities are in place.

This symbol indicates a caution or warning, please refer to the manual.

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## General specification

### Prism and dish

Dish material	316 stainless steel
Prism seal	Silicone rubber
Prism material	Optical glass
Sample surface	8 mm diameter

### Housing

Material	Acrylonitrile Butadiene Styrene
IP rating	IP65 water resistant
Relative humidity	95% RH

### Measurement accuracy

Brix / Refractive Index	±0.2 / ±0.0003
-------------------------	----------------

### Physical

Length	115 mm
Width	54 mm
Height	30 mm
Weight	85g (without batteries fitted)

### Temperature

Storage	-10 to 60°C
Operating	5 to 40°C
Brix measurement	5 to 60°C
Other scale measurement	5 to 40°C

## Automatic temperature compensation

Automatic temperature compensation will correct readings of water and sucrose solutions to 20°C. It conforms to the published ICUMSA 1978 correction tables which covers the ranges 10 to 40°C and 0 to 80° Brix and has been extended to cover 5 to 70°C by using additional data. Although the correction is specifically applicable to pure sucrose solutions, it is also valid for many sugar based food products. Application (scale) specific temperature compensation has been applied to non-sugar related scales using published or experimental data as appropriate.

## Declaration of conformity

According to ISO/IEC 17050-1 & 2 : 2004 Bellingham + Stanley declares that the OPTi Digital Handheld Refractometer (all models) conforms to the following technical requirements:

EMC	Emissions	EN 61326-1:2006 & 2013	CISPR 11:2003, Class B
	AS/NZS CISPR 11		CISPR 11:2003, Class B
	FCC/CFR 47:Part 15		ANSI C63.4:2003, Class B
	Canadian Standard ICES-003:Issue 4		CISPR 22:1997 inc A2:2003
	Immunity	EN 61326-1:2006 & 2013	IEC 61000-4-2:1995 inc A2:2001
		EN 61326-1:2006 & 2013	IEC 61000-4-3:2002 & 2006
Supplementary	The product herewith complies with the requirements of the EMC Directive 2014/30/EU		

## Warranty and customer care

This refractometer is warranted for 12 months after the date of purchase against any manufacturer defect in materials or workmanship. As this refractometer is a precision optical instrument, care must be taken to ensure correct storage, handling and use of the instrument. Failure to do so could invalidate the instrument's warranty.